
Operating Instructions 使用说明书

Recirculating Chiller

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CH800

Contents

Warning Signs

1. Safety Instructions	
2. Application Scenarios	3
3. Unpacking	3
4. Device Setup	5
5. Preparations	5
6. Operating Panel and Display	7
7. Menu Settings	12
8. Temperature Curve Screen	14
9. Operation Settings	14
10. Maintenance and Cleaning	15
11. Error codes	15
12. Technical Specifications	
13. Warranty	

Please read the operating instructions in full before starting up and follow the safety instructions.
The appearance and specifications are subject to change without notice.

WARNING SIGNS

DANGER

Indicates an extremely hazardous situation, which, if not avoided, will result in death, serious injury.

WARNING

Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

CAUTION

Indicates a potentially hazardous situation, which, if not avoided, can result in injury.

NOTICE

Indicates practices which, if not avoided, can result in equipment damage.

1. SAFETY INSTRUCTIONS

General Information:

- Please read the operating instructions carefully before use and follow the safety instructions.
- Keep the operating instructions in a place where it can be accessed by users.
- Ensure that only trained staff can work with the device.
- Follow the safety instructions, guidelines, occupational health, safety and accident prevention regulations.
- Set up the device in a spacious area on an even, stable, clean non-slip, dry and fireproof surface.
- Do not operate the device in explosive atmospheres, or with hazardous substances.
- Protect the device and accessories from bumping and impacting.
- Check the device and accessories for damage before each use. Do not use damaged components.
- Safe operation is only guaranteed with the accessories described in the "Accessories" section.
- The device must only be operated with the original power cable.
- The socket for the power cable must be easily accessible.
- Socket must be grounded (protective ground contact).
- The voltage stated on device label must match the supply voltage.
- The device can only be disconnected from the power supply by pulling out the power plug or the connector plug.
- Disconnect the power plug before attaching or changing any accessories.
- Disconnect the power plug before cleaning, maintenance and transportation of the unit.
- The device must only be opened by trained specialists. The device must be unplugged from the power supply before opening. Some parts inside the device may still cause electrical discharge for some time after unplugging from the power supply.

Coverings or parts that are removable must be put back on the device again to ensure safe operation, for example to keep foreign objects and liquids, etc. from getting into the device.

- The device may only be used as prescribed and as described in these operating instructions. This includes operation by trained personnel.
- When using critical or hazardous materials in your processes, it is recommended to use additional appropriate measures to ensure safety in the experiment. For example, users can implement comprehensive monitoring equipment.
- Process pathogenic material only in closed vessels under a suitable fume hood. Please contact supplier for technical support if you have any question.

If the power switch is not within reach when device is operating, an EMERGENCY STOP switch that can be easily accessed must be installed in the work area.

- The circulating chiller refrigerates and circulates fluid according to specified parameters. There are potential hazards due to high and low temperatures and general electrical hazards while using the device. The user safety cannot be ensured only with design requirements of the device. Other hazards may arise, e.g., by the breakage of the container and reaction with the carrier fluid. It is not possible to consider all eventualities. They remain largely subject to the judgment and responsibility of the operator. For this reason, it become necessary for user to take other precautionary safety measures.
- Use suitable hoses for connection.
- Secure hoses and tubes against slippage and avoid kinks.
- Check hoses, tubes and bath at regular intervals for possible material fatigue (cracks/leaks).
- When device is used for external circulation, extra precaution must be taken for cold fluid leakage due to damaged hose.

**DANGER**

Do not start up the device if:

- It is damaged or leaking
- Cable (not only supply cable) is damaged.
- After a power failure during operation, the device may recover the previous operating mode.
- Transport the device with care.
- Do not transport or empty the bath while it is still cold. Check the temperature of the fluid when emptying the bath.
- Always empty the bath before moving the device.

**NOTICE**

Always empty the bath if the unit is unused for long time.

Fluids:**CAUTION**

Only use the fluids, which fulfill the requirements for safety, health and device compatibility. Be aware of the chemical hazards that may be associated with the bath fluid used. Observe all safety warning for the fluids.

- Depending on the bath fluid used and the type of operation, toxic or combustible vapors can arise. Ensure suitable ventilation.
- Do not use any fluid which may cause dangerous reactions during processing.
- Only use recommended bath fluid. Only use non-acid and non-corroding fluid.

**CAUTION**

Never operate the device without sufficient fluid! Check the fluid level detection at a regular basis.

- Continuous monitoring of the filling level of the bath fluid is required.
- To ensure a sufficient fluid circulation, the bath fluid should be loaded to exceed the yellow line of the filling level indicator. (Recommended fluid volume is 4 to 4.5L)
- Untreated tap water is not recommended. It is recommended to

use distilled water or high purity water (ion ex-changers) and add 0.1g soda (sodium carbonate Na_2CO_3) /liter, to reduce corrosive properties.

Don't use following fluids

- Untreated tap water
- Acids or bases
- Solutions with halides: chlorides, fluorides, bromides, iodides or sulfur
- Bleach (Sodium Hypochlorite)
- Solution with chromates or chromium salts
- Glycerine
- Ferrous water

2.Application Scenarios

• Application :

Used for refrigerating and circulating fluids.

Intended Use: Benchtop device.

Range of use (Indoor environments only) :


- Laboratory
 - College
 - Pharmaceutical
 - Medical/Research institutions
- ### • The device can be used in any areas except:
- Residential areas
- ### • The safety of the user cannot be guaranteed:
- If the device is operated with accessories that are not supplied or recommended by supplier.
 - If the device is operated improperly or in contrary to the supplier's specifications.
 - If the device or the printed circuit board are modified by third parties.

3. Unpacking

- Unpacking :
 - Unpack the device carefully.
 - Any damage should be notified immediately and send damaged note to the shipping agent.
- Packing list :

Item	QTY
CH800 Chiller (Recirculating Chiller)	1
Power cable	1
Inlet & Outlet hose with insulated pipe (inner diameter $\phi 10\text{mm} \times 1.5\text{m}$)	2
Overflow hose (inner diameter $\phi 10\text{mm} \times 1.0\text{m}$)	1
Drain hose (inner diameter $\phi 10\text{mm} \times 1.0\text{m}$) + 1 stopper (outer diameter $\phi 10\text{mm}$)	1
Plastic pipe clamp	6
User manual	1

4. Device Setup

1	Filling lid	10	RS485 interface
2	Handle	11	Power socket
3	Touch screen	12	External Temperature sensor interface
4	Power switch	13	Air outlet
5	Venting grid cover	14	Overflow
6	Filling level indicator	15	Pump connector (OUT)
7	Venting grid	16	Pump connector (IN)
8	Legs	17	Drain
9	USB interface		

5. Preparations

• Setup


- Place the unit on an even, stable, clean, nonslip, dry and fireproof surface.
- Keep at least 20 cm of open space on the front and rear side.
- The place for installation should be large enough and provide sufficient air ventilation to ensure the room does not warm up excessively because of the heat from device radiates to the environment.
- Do not set up the device in the immediate vicinity of heat sources and do not expose to sun light.
- The chiller motor and electronics produce intrinsic heat that is dissipated via the venting grid! Never cover the venting grid!
- Connect the tubing
- Remove the stoppers from the pump connection IN and OUT, overflow and drain.
- Connect hoses to the pump connection IN and OUT and secure the hoses with suitable clamps. The other end of the hose should be connected to the device that needs to be cooled.

- Connect a hose to the “Overflow” port and secure it with clamp. Overflowing fluid should be directed into a vessel via the hose.
- Connect a hose to the “Drain” port and secure it with clamp. Draining fluid should be directed into a vessel via the hose.

NOTICE

When the external system is not necessary, please seal the pump connection IN and OUT with silicone stoppers.

Filling and draining :

- Filling: Open the filling opening lid to fill 4 to 4.5L fluid to the bath and put the lid back.
- Draining: To drain the fluid from the bath, remove the stopper from the drain hose and put the hose in a vessel.
- Connect a hose to the “Overflow” port to discharge excessive fluid to a vessel. The vessel should be placed lower than the overflow port.

NOTICE

Tap water may be unsuitable for operation because the calcium carbonate content may cause calcification. High purity water (from ion exchangers) and distilled or bi-distilled water are unsuitable for operation due to

corrosive properties of these media. High purity water and distillates are suitable as a medium after adding 0.1g soda (Na_2CO_3 , sodium carbonate) per liter of water.

NOTICE

The temperature limit values are adjustable in accordance with the refrigerant used.

NOTICE

Accidental power failure during operation, please wait for 5 minutes before restart.

NOTICE


Preparing suitable bath fluids is required when the working temperature is set under 4°C.

Bath Fluid	Proportion	Minimum Temp. (freezing point)	Recommended Temp. Settings	Volume
Ethylene Glycol (Antifrogen® N)*	5:5 (fluid: water)	-32~-35°C	-10~-20°C	4-4.5L (mixture)
	4:6 (fluid: water)	-22~-25°C	-5~-10°C	4-4.5L (mixture)
	3:7 (fluid: water)	-12~-15°C	0~-5°C	4-4.5L (mixture)
	2:8 (fluid: water)	-2~-5°C	5~0°C	4-4.5L (mixture)
Water	/	4~0°C	Room Temp. ~ 5°C	4-4.5L

*Note: Recommended ratio for Antifrogen® N at different temperature settings

Ethylene glycol solution is not suitable for direct operation and needs to be mixed with distilled water or deionized water in a proper scale.

6. Operating Panel and Display

No.	Designation	Function
A	Bath temperature	Read the actual bath temperature
B	Time	Read the current time
C	"Menu" icon	Press it to display menu option
D	"Temperature curve" icon	Press it to display the bath temperature curve
E	Compressor	"Blue": working; "Gray": not working
F	Cooling fan	"Blue": working; "Gray": not working
G	Pump	"Blue": working; "Gray": not working
H	Temperature sensor	"Blue": external temperature sensor "Gray": internal temperature sensor
I	"Set" icon	Press it to display the operation settings options
J	"On/Off" icon	Switch on/off the chiller
K	"Timer Setpoint"	Press it to set the time: 0-30h
L	Timer display	Displays the remaining working time
M	"Temp Setpoint"	Set the temperature: -20°C~ room temperature

Before starting up, make sure the voltage specified on the device label matches the available supply voltage.

The power socket used must be grounded (protective earthing).

If above conditions are met, the device is ready for operation after connecting the power.

If these conditions are not met, safe operation is not guaranteed, and the device could be damaged.


Observe the ambient conditions (temperature, humidity, etc.) listed under “Technical Parameters”.

Power on

1. Switch on the power switch.
2. Then the device enters standby status and is ready for operation.
3. Press “Temp Setpoint” to set the temperature.
4. Press “Timer Setpoint” to set the time.
5. Make sure all settings are correct, press “On/Off” icon to start the device.

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
[Cancel](#) [Save](#)

- **Temperature Unit and Screen Brightness**

Temp.Unit

Temperature unit ☒ °C ☐ °F

Screen brightness 50 100

[Back](#)

• Temperature Calibration

1. In the home screen, set the temperature value at 0°C. When it stabilizes to 0°C, use an external temperature probe to measure the actual temperature of the fluid, and input the measured temperature value to the “Calibrated Temp.” box.
2. In the home screen, set the temperature value at -15°C. When it stabilizes to -15°C, use an external temperature probe to measure the actual temperature of the fluid, and input the measured temperature value to the “Calibrated Temp.” box.
3. Then click “Calibrate” to perform the settings.

Calibration

Setpoint	Set Temp.	Current Temp.	Calibrated Temp.	Confirm
Point 1	0	0	0	Yes
Point 2	-15	-15	-15	Yes

• Device Information

Device Information	
Software version	Ver2.0
Firmware version	Ver2.0
Model Number	CH800

Back


• RS485 Interface

RS485 Interface	
Baud rate	115200 ▼
RS485 address	0001
Check bit	None ▼
Data bit	8 ▼
Stop bit	1 ▼

Cancel Save

8. Temperature Curve Screen

- Temperature curve :

No.	Designation	Remarks
A	Temp. value	Read the temperature
B	Time value	Read the time
C	"Back" button	Back to the Home Screen

9. Operation Settings

• Operation Settings:

Operation Settings

A

Sensor

☒ Internal ☐ External

B

Pump

☒ ON ☐ OFF

C

Timer

☒ ON ☐ OFF

D

Back

No.	Designation	Remarks
A	"Sensor" setting	Select Internal/External temperature sensor
B	"Pump" setting	Switch on/off the pump
C	"Timer" setting	Switch on/off to set the time
D	"Back" button	Back to the home screen

10. Maintenance and Cleaning

To avoid contamination, it is essential to change the bath fluid regularly. If water is used for bath fluid, we recommend adding water bath protective agent. The protective agent stops the growth of algae, bacteria and other microorganisms to ensure the bath fluid can be used for long time. To keep the full cooling performance, the dust filter of the cooler must be checked regularly and cleaned if necessary.

- Switch off the device and disconnect power cable.
- Detach the venting grid cover.
- Take out the dust filter.
- Clean the dust filter with a vacuum cleaner or wash it in the water and dry it before assembly.

NOTICE

Don't clean the compressor surface with sharp/hard objects.

Cleaning

- Disconnect the power prior to cleaning.
- Use only washing agents which have been approved by supplier to clean the device: water containing surfactant or isopropyl alcohol.
- Wear protective gloves while cleaning the device.
- Electrical devices may not be placed in the cleansing agent while cleaning.
- Do not allow moisture to get into the device when cleaning.
- If use a different cleaning or decontamination method that is not recommended by the supplier, the user must check with the supplier that this method does not damage the device.

Repair

- Please send in device for repair only after it has been cleaned and is free from any materials which may pose a health hazard.
- For repair, please return the device in its original packaging. If the original packaging is missing, please use suitable packaging for transportation.

11.Error Codes

Any malfunctions during operation will be identified by an error message on the display.

Proceed as follows in such cases:

1. Switch off the device
2. Restart the device

Error message	Probable Cause	Suggested Action
Calibration error	Entered the wrong calibration temperature value	Check the parameters and re-enter
Internal temperature sensor fault	Damaged sensor or loose wires	Contact customer service department
External temperature sensor fault	Damaged sensor or loose wires	Check if the external temperature sensor is plugged properly into the device.
Fluid temperature is too high	Fluid temperature is higher than 40°C	Switch the device off and let the fluid cool down
Over current	Damaged compressor	Power off the device, wait for 10 minutes, and restart the device.

If the actions described above fails to resolve the fault or another error message is displayed, then contact customer service department.

12. Technical Parameters

Operating voltage		VAC	AC100~120V/AC200~240V
Frequency		Hz	50Hz/60Hz
Max. input power		W	750
Working temperature range		°C	-20°C ~ room temperature (25°C)
Temperature stability		K	±0.3
Temperature control		°C	PID
Temperature measurement, absolute accuracy		K	±0.5
Temperature setting		/	Touch screen
Temperature setting resolution		K	0.1
Temperature display		/	Touch screen
Temperature display resolution		K	0.1
Cooling capacity	+20°C	W	750
	+10°C		650
	0°C		500
	-10°C		400
	-20°C		200
Refrigerant		/	R290
Max. Pump pressure		bar	1.0
Max. pump flow rate		l/min	22
Bath volume		L	5.0L(Max. 5.3L)
Fluid maximum viscosity		mm2/s	50
Interface		/	USB; RS485
Permitted ambient temperature		°C	+5~+35°C
Permitted ambient humidity		%	80
Dimension (W x D x H)		mm	235*468*535mm
Weight		kg	35

Refrigerant must be disposed of in accordance with local and national regulations. The technical parameter is subject to change without notice!

13.Warranty

We guarantee that our scientific instruments adhere to the most rigorous engineering and quality standards. This instrument is warranted to be free from defects in materials and workmanship under normal use and service, for a period of 24 months from the date of dispatch. The warranty is extended only to the original purchaser. For claims under the warranty, please contact your local supplier. After the warranty period expires, the manufacturer retains the right to invoice the cost price for the repair or maintenance of a faulty device, along with any associated service fees.

Scope of Warranty

The following conditions are not covered under the warranty.

- Faults or damage caused by negligence, improper installation, improper operation, or failure to use and maintain the machine in accordance with the instructions in this operating manual.
- Issues caused by unauthorized disassembly or modification.

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