# Operating Instructions 使用说明书

LED Digital Rocking Shaker  
LED数控翘板摇床  
Blots / SK308

# Contents

Preface	2
Service	2
Warranty	2
1. Safety Instructions	3
2. Proper Use	4
3. Inspection	4
3.1 Receiving Inspection	4
3.2 Listing of Items	5
4. Control and Display	5
4.1 Control	5
4.2 Display	6
5. Operation	7
6. Faults	8
7. Maintenance and Cleaning	8
8. Associated standards and regulations	9
9. Technical Data	10
10. Accessories	11

## **Preface**

Welcome to the "Operating Instructions of Rocking Shaker". Users should read this manual carefully, follow the instructions and procedures, and beware of all the cautions when using this instrument.

## **Service**

When help needed , you can always contact us for technical support. Please provide the customer care representative with the following information:

- Series number ( on the rear panel )
- Description of problem ( i.e., hardware or software )
- Your contact information

## **Warranty**

You have purchased an original laboratory machine which meets the highest engineering and quality standards. In accordance with our warranty conditions, the warranty period is 24 months from the date of invoice. It shall not apply to any product or parts which have been damaged on account of improper installation, improper connections, misuse, accident or abnormal conditions of operation.

For claims under the warranty please contact your local dealer. You may also send the machine directly to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

# 1. Safety Instructions

## Warning!

- Read the operating instructions carefully before use.
- Ensure that only trained staff works with the instrument

## Protective ground contact !

- Ensure the socket must be grounded (protective ground contact) before use.
- When working wear personal safety guards to avoid the risk of:
  - Splashing liquids
  - Broken glass containers
  - Body, hair, clothes and jewelry are involved
- Follow the safety instructions, guidelines and accident prevention regulations.
- Do not touch the running parts, prevent the hair and clothes involved. When moving instrument care not rolling your fingers.
- Set up the instrument on a stable, clean, non-slip, dry and fireproof surface. Do not operate the instrument in explosive atmospheres, with hazardous substances or under water.
- When setting speed, please pay attention to the vessels on shaking table to avoid the sample splashing.  
If the instrument does not run smoothly, please decrease the motor speed.
- Firmly secure the accessories and vessels in place to avoid damage or risk.

- Processing flammable samples will be dangerous. Only process samples that will not react dangerously.
- Use the standard accessories listed in the “accessories” section, and follow the instructions to use accessories to ensure safety. Please switch off the power before assembly of accessories, confirm the instrument and accessories are intact before switching on each time.
- Only professionals can open the instrument. Please pull out the power plug before opening the instrument.
- The voltage stated on the nameplate must correspond to the mains voltage.
- Do not cover the instrument during running.
- Keep the instrument away from the high magnetic field.

## **2. Proper Use**

The instrument is designed for shaking and mixing liquid in schools, laboratories or factories. It can be installed with a variety of accessories and suitable for kinds of tubes. This instrument is not suitable for using in residential areas or other constraints mentioned in Chapter “ Safety Instructions ” .

Do not use the accessories recommended by the manufacturer, or failure to use the instructions, may be caused unsafe situation.

## 3. Inspection

### 3.1 Receiving Inspection

Unpack the instrument carefully and check for any damages which may happen during transport. If any damage happens, please contact us for technical support.

**Note:**


If there is any apparent damage to the system, please do not connect to the power line.

### 3.2 Listing of Items

<b>Item</b>	<b>Qty.</b>
Main unit	1
Power adapter	1
Non-slip mat	1
Fixed string	6
Operating instruction	1


## 4. Control and Display

### 4.1 Control

Items	Descriptions
Speed control knob	Set target speed. Push on/off the knob to start/stop the shaking function.
Timer control knob	Set target time. Push on/off the knob to start/stop the shaking function.
Speed display	Display the set and real speed values.
Time display	Display the set and remain time values.
Power switch	Switch ON/OFF instrument.

## 4.2 Display

Display	Descriptions
Er 1	Means motor blocked .
Er 2	Means motor overspeed.
Display area	When the value on display area is flashing, the current display value is setting value; When value stops flashing, the current display value is running value.

## 5. Operation

- Place the instrument on safe and stable surface, ensure the required operating voltage and power supply voltage matched, and the socket must be earthed reliably.
- Connect the power cable.
- Switch "ON" of the instrument.
- The instrument begins self-checking.
- When initialization is over, LED displays the last running values of speed and time.
- Rotate the speed control knob to set target speed.
- Press the speed control knob to start shaking function without timer. Press the speed control knob again to stop shaking function.
- Press the timer control button to start shaking, and rotate the timer control knob to set target time. Before the countdown time is over, the shaking can be stopped at any time by pressing the speed control knob.

When the timer is closed by pressing the timer control knob, press the timer control knob again, the shaking will start again and the timer will restart countdown. When the timer reaches zero, the unit will be automatically halted and an alert will sound.

If these above operations are normal, the instrument is ready to operate. If not, the instrument may be damaged during transportation, please contact technical support of the manufacturer/supplier.

**Note:**

1. The speed should be adjusted by speed control knob slowly to make the instrument run smoothly.
2. The speed value can be altered during shaking. The instrument will start to shake with new setting speed value.

## **6. Faults**

- When switched ON, the instrument does not work.  
Check whether the power cable is connected.  
Check whether the fuse is broken.
- The running speed can't reach to the set value.  
Please check whether the load of the shaker exceeds the max. load capacity.

## **7. Maintenance and Cleaning**

- Proper maintenance can keep instrument working properly and lengthen its lifetime.
- Do not spray cleanser into the instrument when cleaning.

- Do not remove the power line when cleaning. Only use recommended cleansers:

Dyes	Isopropyl alcohol
Construction materials	Water containing tenside / Isopropyl alcohol
Cosmetics	Water containing tenside / Isopropyl alcohol
Foodstuffs	Water containing tenside
Fuels	Water containing tenside

- For cleansers not listed, please consult the manufacturer's after-sales service department first.
- Before using other method for cleaning or decontamination, the user must ascertain with the manufacturer that this method will not damage the instrument. Wear the proper protective gloves during cleaning of the instrument.

**Note:**

- Electronic device cannot be clean with cleanser.
- If maintenance service required, the instrument must be cleaned first to avoid pollution of hazardous substances.
- If the instrument will not use for a long time, please switch it off and place it in a dry, cleans,room temperature and stable location.

## 8. Associated standards and regulations

Construction in accordance with the following safety standards:

EN 61010-1

UL 3101-1

CAN/CSA C22.2(1010-1)

EN 61010-2-10

Construction in accordance with the following EMC standards:

EN 61326-1

## 9. Technical Data

Items	Specifications
	Rocking Shaker SK308
Voltage [VAC]	100-240
Frequency [Hz]	50/60
Power [W]	40
Tile Angel [°]	9°
Max. Load capacity(with platform) [kg]	10
Motor Type	Brushless DC motor
Speed range [rpm]	10-70
Speed display	LED
Timer display	LED
Timer range [min]	1min-99h59min
Run type	Continuous / timing operation
Overall Dimensions [mm]	452x407x175
Weight [kg]	9.7
Permissible ambient temperature [°C]	5-40
Permissible relative humidity	80%
Protection class acc. to DIN EN60529	IP21

## 10. Accessories

	Cat. No.	Descriptions
	SK308-2-10	Dish platform with non-slip mat and fixed string(6 pcs)
	MI0103002-101	Universal attachment
	MI0103002-102-1	Fixing clip attachment
	MI0103002-102a	Fixing clip for flask volume 25 ml
	MI0103002-102b	Fixing clip for flask volume 50 ml
	MI0103002-102c	Fixing clip for flask volume 100 ml
	MI0103002-102d	Fixing Clip for flask volume 250ml
	MI0103002-102e	Fixing Clip for flask volume 500ml

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