



NEW

American Dilation System with CleanGuide™

Reusable Polyvinyl Dilators

This over-the-wire system of radiopaque dilators provides a full range of diameter sizes for the dilatation of esophageal strictures while the CleanGuide™ guidewire is designed to improve infection prevention by eliminating potential cross contamination that could occur with reusable guidewires.

American Dilator

NEW CleanGuide™ Disposable
Marked Spring Tip Guidewire

Features and Benefits

Effective Dilation

The Reusable Polyvinyl Dilators of the American Dilation System have features designed to provide dilation that is both clinically and cost-effective

- Graduated sizes
- Radiopaque dilators
- Calibrated American and European measurement system
- Short 4cm distal tip
- 80cm overall length

Infection Prevention

The sterile CleanGuide™ Marked Spring Tip Guidewire helps reduce the risk of potential cross contamination that could occur from using a reusable guidewire.*

Patient Safety

Over time, the overuse of reusable guidewires can eventually cause the tip to bend, fatigue, discolor, become loose or crack which may cause harm. The single-use, disposable design of CleanGuide™ guidewires helps mitigate this concern.

Facility Safety

CleanGuide™ is a single-use product that eliminates the need to reprocess the wire. By eliminating the reprocessing we limit the exposure of hazardous material to the staffs cleaning the wire and limit the ability for the wire to spring open during reprocessing since it is a stiff wire.

*Looking for a dilation
alternative?*

**CONMED also offers
Eliminator® PET
Balloon Dilators.**

*The American Dilation System is also available with a reusable marked spring tip guidewire. CleanGuide™ is packaged sterile to offer an additional level against potential cross contamination.

To take the next step to safe and effective dilation, call **1-866-4CONMED** to arrange a demonstration with your local CONMED Sales Representative, or visit **CONMED.com** for more information.

Each American Dilation System* Includes:

15 ea.	Tapered over-the-wire dilators	1 ea.	Cleaning brush
1 ea.	Marked spring tip guidewire	1 ea.	Carrying case
2 ea.	Cleaning adaptor		



Ordering Information

American Dilation System

Non-Sterile, Reusable

Description	Size	Unit	Catalog Number
American Dilatation System		1/ box	000250*
Tapered Over-The-Wire Dilator	15 F	1/box	000270
Tapered Over-The-Wire Dilator	18 F	1/box	000288
Tapered Over-The-Wire Dilator	21 F	1/box	000271
Tapered Over-The-Wire Dilator	24 F	1/box	000289
Tapered Over-The-Wire Dilator	27 F	1/box	000272
Tapered Over-The-Wire Dilator	30 F	1/box	000273
Tapered Over-The-Wire Dilator	33 F	1/box	000274
Tapered Over-The-Wire Dilator	36 F	1/box	000275
Tapered Over-The-Wire Dilator	39 F	1/box	000290
Tapered Over-The-Wire Dilator	42 F	1/box	000277
Tapered Over-The-Wire Dilator	45 F	1/box	000278
Tapered Over-The-Wire Dilator	48 F	1/box	000279
Tapered Over-The-Wire Dilator	51 F	1/box	000280
Tapered Over-The-Wire Dilator	54 F	1/box	000281
Tapered Over-The-Wire Dilator	60 F	1/box	000291

Accessories for American Dilators

NEW	CleanGuide™ Disposable Marked Spring Tip Guidewire*	10/box	DIS150
	American Dilator Carrying Case	1/box	000249
	American Dilator Cleaning Brush (single use)	3/box	000336
	Cleaning Adaptor (single use)	2/box	000337

*The American Dilation System is also available with a reusable marked spring tip guidewire.

CONMED Corporation
525 French Road, Utica, NY 13502

International: 727-214-3000
Toll Free: 1-866-4CONMED

CONMED.com
customer_service@conmed.com