

24k[®]

High Performance Brilliantly Controlled

You make the choice... the 24k[®] delivers with superior performance



HIGH PERFORMANCE DEAM HIGH PERFORMANCE D

The 24k[®] delivers high performance with custom features that allow you to get the job done.

- › Independent pressure and flow controls allow the 24k[®] Pump to keep the joint clear of debris even at lower pressure settings.
- › Independent rotor control allows the 24k Pump to operate as a simple inflow-only pump or a complete arthroscopic fluid management system.
- › Four Presets – Knee, Shoulder, Small Joint and Hip simplify setting up the pump.
- › Custom Wave Mode in the Knee Preset for excellent joint clearing even with a double stopcock scope sheath.
- › Simple intuitive tubing set up. Hang the cassettes, spike the bags, prime and go.
- › Independent shaver sensor allows the 24k Pump to control suction with any ConMed Linvatec or competitive shaver handpiece.
- › Fully autoclavable multifunction remote control gives you complete control of the pump.



DEMANDED... DELIVERED.

functions and flexible

› Higher Suction › Higher Outflow
› Wave Mode › Shaver Flow Compensator
› TIPS Performance › Clearfield Performance



24k® PUMP ORDERING INFORMATION

Cat. No.	Description
24k	24k® Pump
24k900	24k Shaver Sensor
24k191	24k Shaver Sensor Extension Cable
C8115	24k Remote
24k100	24k Outflow/Shaver Suction Tube Set
24kTIPS	24k TIPS Tube Set
10k100	10k® Arthroscopy Tube Set <i>(Inflow Only)</i>
10k100D	10k Arthroscopy Tube Set for Dyonics
10k150	10k Arthroscopy Tube Set with Suction Line
10k500	10k Laparoscopic Tubing Set
10k505	10k Laparoscopic Tube Set with Probe
10k600	Eco-Flow™ Day-use Pump Tubing <i>(Must be used with 10k605)</i>
10k605	Eco-Flow™ Single-use Patient Tubing <i>(Must be used with 10k600)</i>
C7104	AC Power Cord 115VAC
C7105	AC Power Cord 230 VAC

ConMed Linvatec Fluid Management System

11311 Concept Boulevard
Largo, Florida 33773-4908

Local 727-392-6464
Toll Free 800-237-0169

Website www.linvatec.com
Email customer_service@linvatec.com